

MICRODENT CAD-CAM

PRODUCT CATALOGUE

www.microdentsystem.com



CONTENTS

04

Microdent CAD-CAM

- 05 Introduction
- 06 Intraoral scanbody
- 07 Laboratory scanbody
- 08 Digital replica
- 09 Straight Dome Ti Base
- 10 Dynamic Dome Ti Base
- 11 Screw
- 12 Straight and angled Rotating Mini Capitel Transepithelial Abutment
- 13 Mini Hexagonal Capitel Transepithelial Abutment

28

Universal external connection implant library

- 29 Direct implant procedure
- 31 Dome Ti Base Procedure
- 34 Direct Mini Hexagonal Capitel Abutment Procedure
- 36 Ti Base Mini Hexagonal Capitel Abutment Procedure
- 38 Direct Rotating Mini Hexagonal Capitel Abutment Procedure
- 40 Ti Base Rotating Mini Capitel Abutment Procedure

57

Ektos hexagonal internal connection implant library

- 58 Direct implant procedure
- 60 Dome Ti Base Procedure
- 64 Direct Mini Hexagonal Capitel Abutment Procedure
- 66 Ti Base Mini Hexagonal Capitel Abutment Procedure
- 68 Direct Rotating Mini Hexagonal Capitel Abutment Procedure
- 70 Ti Base Rotating Mini Capitel Abutment Procedure

14

Microdent external connection implant library

- 15 Direct implant procedure
- 17 Dome Ti Base Procedure
- 20 Direct Mini Hexagonal Capitel Abutment Procedure
- 22 Ti Base Mini Hexagonal Capitel Abutment Procedure
- 24 Direct Rotating Mini Hexagonal Capitel Abutment Procedure
- 26 Ti Base Rotating Mini Capitel Abutment Procedure

42

Genius conical internal connection implant library

- 43 Direct implant procedure
- 45 Dome Ti Base Procedure
- 49 Direct Mini Hexagonal Capitel Abutment Procedure
- 51 Ti Base Mini Hexagonal Capitel Abutment Procedure
- 53 Direct Rotating Mini Hexagonal Capitel Abutment Procedure
- 55 Ti Base Rotating Mini Capitel Abutment Procedure

72

Summary of galleries and procedures

MICRODENT CAD-CAM






MICRODENT CAD-CAM

INTRODUCTION

Introducing the **MICRODENT CAD-CAM** system. In this catalogue you will find all the references corresponding to the prosthetic and implant attachments included in our libraries: Direct connection, Dome Ti base or Mini Capitel Transepithelial Abutment

At **Microdent**, we are convinced that the greatest benefit that a prosthesis can bring to the durability of implants in the mouth is its passivity. For this reason, we use cement-screw-retained prostheses; cemented on our original attachments and screwed with maximum passivity. With this technique, we also ensure that we always work with the same material, titanium, both in the implant and the connection to the implant as well as the screw that secures the assembly.



INTRAORAL SCANBODY

Microdent's new intraoral scanbody is designed for digital recording of the implant position by means of an intraoral scanner in the clinic. With a superb geometric design and high surface quality, it provides high-precision scanning results. Manufactured in aluminium that prevents damage to the implant connection without losing fit or precision, and with a sandblasted surface that allows easy and quick reading, the exact location of the implant position avoids complicated specific pre-treatments. The radiological opacity makes it possible to check the scanbody adjustment to the implant avoiding errors from the beginning of the work. Two models of scanbody, short and long, make it easy to adapt to the needs of different patients.

Microdent's intraoral scanbody can be used for direct restorations of implants with angled and straight prostheses, indirect restorations of implants via straight interfaces (Dome Ti base) or angled interfaces (dynamic Dome Ti base) and overdentures. Always with the passive fit of our prosthesis as the main objective.




LABORATORY SCANBODY

Microdent's new laboratory scanbody is designed for digital recording of the implant position by means of an laboratory scanner in the laboratory. With a superb geometrical design and high surface quality, it provides high-precision scanning results thanks to the fit over the screwed-in intermediate analogues (replicas). Anti-rotational, made of aluminium and featuring a sandblasted

surface that allows easy and quick reading as well as the exact location of the implant position.


The tendency to reduce the size of the scanbodies in order to achieve more anatomically compatible morphologies must never compromise the accuracy of the scan, and priority must always be given to the level of precision of the digital technology.

Microdent's laboratory scanbody can be used for direct restorations of implants with angled and straight prostheses, indirect restorations of implants via straight interfaces (Dome Ti base) or angled interfaces (dynamic Dome Ti base) and overdentures. Always with the passive fit of our prosthesis as the main objective.



DIGITAL REPLICA (digital analogue)

Microdent's digital replicas are manufactured in grade V titanium in order to exactly replicate the geometry of the implant connection in a 3D manufactured model, facilitating its morphology and excellent positioning. It is essential that the positioning is accurate to avoid errors at the time of making the prosthesis. For this purpose, the analogues have a retentive zone at the apex which is secured by a screw that holds the replica in the model avoiding vertical movements, and three flat faces that allow lateralisations to be controlled. They are anodised in the colour of the implant system for easy identification.




STRAIGHT DOME TI BASE

Microdent's straight Dome (Ti base) interfaces are made of grade V titanium and provide a bond in customised digital restorations. The titanium base is bonded to the prosthesis for a perfect connection to the implant or Mini Capitel transepithelial abutment, depending on the case. Available in short and long lengths, the latter can be cut, with different gingival heights depending on the system. Rotational, with a slightly conical prosthetic section to facilitate the insertion of multiple restorations, and anti-rotational for individual restorations, the Dome interface can be adapted to the varying characteristics of the patients.

They are sandblasted for better bonding/retention and gold anodised for better translucency.

Dome abutments make it possible to create the ideal prosthesis, combining the advantages of a single-unit prosthesis with a


bridge made of zirconium, PMMA (polymethylmethacrylate) or other materials.




DYNAMIC DOME TI BASE

Microdent's angled Dome (Ti base) interfaces are made of grade V titanium and provide a bond in customised digital restorations. The titanium base is bonded to the prosthesis for a perfect connection to the implant or Mini Capitel transepithelial abutment, depending on the case. Available in short and long lengths, with different gingival heights depending on the system, rotating and anti-rotational for multiple or individual restorations, the Dome interface can be adapted to the varying characteristics of the patients. They are sandblasted for better bonding/retention and gold anodised for better translucency. The main feature of the dynamic Dome abutments is to provide a solution for angled implants by means of the Cad-Cam, which is able to correct up to 25 degrees of divergence. They are joined together with the Torx screw.

Dome abutments make it possible to create the ideal prosthesis, combining the advantages of a single-unit prosthesis with a bridge made of zirconium, PMMA (polymethylmethacrylate) or other materials.




 See table of angles according to heights.




SCREW

Microdent screws are made of grade V titanium and are supplied in sets with the attachments. All of them have a universal 1.20 mm universal screwdriver connection, except for the Torx model (dynamic Domes) which is made by Microdent.


The Torx screw is only used to fix the dynamic Domes. These elements are recommended for implant restorations in angled positions. The geometry of this connection (Torx) facilitates screwing and unscrewing in these unfavourable positions.



 Maximum torque standard screw: 30Nw. Maximum torque Torx screw: 20Nw


ANGLED AND ROTATING STRAIGHT MINI CAPITEL TRANSEPITHELIAL ABUTMENT

Mini Capitel abutments are Microdent's transepithelial abutments made of grade V titanium and anodised to maximise prosthetic imitation even in thin biotypes, allowing prosthetic solutions to be even more individualised. Mini Capitel transepithelial abutments are designed for multiple screw-retained restorations. The system has different gingival heights and angulations, between 17° and 30°, to correct divergences between implants. They are supplied together with their screw and conveyor and are attachments compatible with other existing implant systems. Their reduced size means they can be used with a complete range of prosthetic applications.



HEXAGONAL STRAIGHT MINI CAPITEL TRANSEPIHELIAL ABUTMENT

Mini Capitel abutments are Microdent's transepithelial abutments made of grade V titanium and anodised to maximise prosthetic imitation even in thin biotypes, allowing prosthetic solutions to be even more individualised. The anti-rotational Mini Capitel transepithelial abutments are designed for individual screw-retained restorations. The system has different gingival heights. They are supplied together with their conveyor and an 1.8 metric screw to render the restoration more robust. Their dimensions allow a full range of prosthetic applications and are ideal for the One Abutment One Time technique.







MICRODENT EXTERNAL CONNECTION IMPLANT LIBRARY




PROCEDURE DIRECT IMPLANT

GALLERY: MICRODENT-1400_HEX_EXT_D





Diameter Platform	Nomenclature	Connection type (Anti-rotational, Rotating or Angled)	Microdent Scanbody	Microdent Replica
 Ø3.50mm	1401 3.5 LAB (Laboratory)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating) E (Engaging - Angled) NE (Non-engaging - Angled)	CMSSBL35	CMSRID35
	1401 3.5 IO C (Short Intraoral)		CMSSBI35C	
	1401 3.5 IO L (Long Intraoral)		CMSSBI35L	
 Ø4.20mm	1402 4.2 LAB (Laboratory)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating) E (Engaging - Angled) NE (Non-engaging - Angled)	CMSBL42	CMRID42
	1402 4.2 IO C (Short Intraoral)		CMSBI42C	
	1402 4.2 IO L (Long Intraoral)		CMSBI42L	
 Ø5.10mm	1403 5.1 LAB (Laboratory)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating) E (Engaging - Angled) NE (Non-engaging - Angled)	CMSBL51	CMRID51
	1403 5.1 IO C (Short Intraoral)		CMSBI51C	
	1403 5.1 IO L (Long Intraoral)		CMSBI51L	
 Ø5.60mm	1404 5.6 LAB (Laboratory)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating) E (Engaging - Angled) NE (Non-engaging - Angled)	CMKSBL56	CMKRID56
	1404 5.6 IO C (Short Intraoral)		CMKSBI56C	
	1404 5.6 IO L (Long Intraoral)		CMKSBI56L	

FLOW DIRECT IMPLANT





PROCEDURE WITH DOME TI BASE


GALLERY: MICRODENT-1400_HEX_EXT_TB

Diameter Platform	Nomenclature	Connection type (Anti-rotational, Rotating or Angled)	Ti Base Straight dome	Ti Base Dynamic dome
 Ø3.50mm	1401 3.5 LAB (Laboratory)	E (Engaging - Anti-rotational)	Short CMSDOM35H	
	1401 3.5 IO C (Short Intraoral)	NE (Non-engaging - Rotating)	Short CMSDOM35R Long CMSDOM35HL Long CMSDOM35RL	Short CMSDOMD35H Short CMSDOMD35R
	1401 3.5 IO L (Long Intraoral)	E (Engaging - Angled) NE (Non-engaging - Angled)		Long CMSDOMD35HL Long CMSDOMD35RL
 Ø4.20mm	1402 4.2 LAB (Laboratory)	E (Engaging - Anti-rotational)	Short CMSDOM42H	
	1402 4.2 IO C (Short Intraoral)	NE (Non-engaging - Rotating)	Short CMSDOM42R Long CMSDOM42HL Long CMSDOM42RL	Short CMSDOMD42H Short CMSDOMD42R
	1402 4.2 IO L (Long Intraoral)	E (Engaging - Angled) NE (Non-engaging - Angled)		Long CMSDOMD42HL Long CMSDOMD42RL
 Ø5.10mm	1403 5.1 LAB (Laboratory)	E (Engaging - Anti-rotational)	Short CMSDOM51H	
	1403 5.1 IO C (Short Intraoral)	NE (Non-engaging - Rotating)	Short CMSDOM51R Long CMSDOM51HL Long CMSDOM51RL	Short CMSDOMD51H Short CMSDOMD51R
	1403 5.1 IO L (Long Intraoral)	E (Engaging - Angled) NE (Non-engaging - Angled)		Long CMSDOMD51HL Long CMSDOMD51RL
 Ø5.60mm	1404 5.6 LAB (Laboratory)	E (Engaging - Anti-rotational)	Short CMSDOM56H	
	1404 5.6 IO C (Short Intraoral)	NE (Non-engaging - Rotating)	Short CMSDOM56R Long CMSDOM56HL Long CMSDOM56RL	Short CMSDOMD56H Short CMSDOMD56R
	1404 5.6 IO L (Long Intraoral)	E (Engaging - Angled) NE (Non-engaging - Angled)		Long CMSDOMD56HL Long CMSDOMD56RL

FLOW WITH STRAIGHT DOME TI BASE





FLOW WITH DYNAMIC DOME TI BASE



PROCEDURE

DIRECT MINI HEXAGONAL CAPITAL ABUTMENT


GALLERY: [MICRODENT-1430_Capitel_E_Mini_D](#)

Diameter Platform	Nomenclature	Connection type (Anti-rotational or Rotating)	Scanbody Mini Capitel	Replica Mini Capitel
 Ø4.20mm	1431 Mini 4X LAB (Laboratory)		CCAPNSBL4XH	
	1431 Mini 4X IO C (Short Intraoral)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating)	CCAPNSBI4XHC	CCAPNRID4XH
	1431 Mini 4X IO L (Long Intraoral)		CCAPNSBI4XHL	
 Ø5.10mm	1432 Mini 5X LAB (Laboratory)		CCAPNSBL51H	
	1432 Mini 5X IO C (Short Intraoral)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating)	CCAPNSBI51HC	CCAPNRID51H
	1432 Mini 5X IO L (Long Intraoral)		CCAPNSBI51HL	



FLOW

DIRECT MINI HEXAGONAL CAPITAL ABUTMENT





Microdent always recommends the use of the intraoral scanbody for a better fit between the prosthesis and the Mini Capitel Abutment using the one abutment one time technique.

PROCEDURE

TI BASE MINI HEXAGONAL CAPITAL ABUTMENT

GALLERY: [MICRODENT-1430_Capitel_E_Mini_TB](#)

Diameter Platform	Nomenclature	Connection type (Anti-rotational or Rotating)	Ti Base Straight dome
 Ø4.20mm	1431 Mini 4X LAB (Laboratory)		
	1431 Mini 4X IO C (Short Intraoral)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating)	CUTSNPC4XH
	1431 Mini 4X IO L (Long Intraoral)		
 Ø5.10mm	1432 Mini 5X LAB (Laboratory)		
	1432 Mini 5X IO C (Short Intraoral)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating)	CUTSNPC51H
	1432 Mini 5X IO L (Long Intraoral)		



FLOW TI BASE MINI HEXAGONAL CAPITAL ABUTMENT

Platform Diameter

■ Ø4.20mm

■ Ø5.10mm



Straight (heights)

- 1mm → CMCAPN4201H
→ CMCAPN5101H
- 2mm → CMCAPN4202H
→ CMCAPN5102H
- 3mm → CMCAPN4203H
→ CMCAPN5103H
- 4mm → CMCAPN4204H
→ CMCAPN5104H



Scanbodies on Mini Capital Transepithelial Abutment

Intraoral scanbody

Laboratory scanbody

Short

- CCAPNSBI4XHC
- CCAPNSBI51HC



- CCAPNSBL4XH
- CCAPNSBL51H



Long

- CCAPNSBI4XHL
- CCAPNSBI51HL



Dome Ti Base for Mini Capital

Straight



- CUTSNPC4XH
- CUTSNPC51H



PROCEDURE


DIRECT ROTATING MINI CAPITEL ABUTMENT

GALLERY: [MICRODENT-1440_Capitel_NE_Mini_D](#)

Diameter Platform	Nomenclature	Connection type (Rotating)	Scanbody Mini Capitel	Replica Mini Capitel
 Ø4.20mm	1441 Mini 4X LAB (Laboratory)	NE (Non-engaging - Rotating)	CCAPNSBL4X	CCAPNRID4X
	1441 Mini 4X IO C (Short Intraoral)		CCAPNSBI4XC	
	1441 Mini 4X IO L (Long Intraoral)		CCAPNSBI4XL	
 Ø5.10mm	1442 Mini 5X LAB (Laboratory)	NE (Non-engaging - Rotating)	CCAPNSBL51	CCAPNRID51
	1442 Mini 5X IO C (Short Intraoral)		CCAPNSBI51	
	1442 Mini 5X IO L (Long Intraoral)		CCAPNSBI51	



FLOW DIRECT ROTATING MINI CAPITEL ABUTMENT





Microdent always recommends the use of the intraoral scanbody for a better fit between the prosthesis and the Mini Capitel Abutment using the one abutment one time technique.

PROCEDURE

TI BASE ROTATING MINI CAPITAL ABUTMENT

GALLERY: [MICRODENT-1440_Capitel_NE_Mini_TB](#)


Diameter Platform	Nomenclature	Connection type (Rotating or Angled)	Ti Base Straight dome	Ti Base Dynamic dome
 Ø4.20mm	1441 Mini 4X LAB (Laboratory)	NE (Non-engaging - Rotating)	CUTSNPC4X	
	1441 Mini 4X IO C (Short Intraoral)			
	1441 Mini 4X IO L (Long Intraoral)	NE (Non-engaging - Angled)		CUTSNDPC4X
 Ø5.10mm	1442 Mini 5X LAB (Laboratory)	NE (Non-engaging - Rotating)	CUTSNPC51	
	1442 Mini 5X IO C (Short Intraoral)			
	1442 Mini 5X IO L (Long Intraoral)	NE (Non-engaging - Angled)		CUTSNDPC51




FLOW TI BASE ROTATING MINI CAPITAL ABUTMENT

Platform Diameter

- Ø4.20mm
- Ø5.10mm



Scanbodies on Mini Capital Transepithelial Abutment



Dome Ti Base for Mini Capital





**UNIVERSAL EXTERNAL
CONNECTION IMPLANT LIBRARY**




PROCEDURE DIRECT IMPLANT

GALLERY: MICRODENT-1410_UNIVERSAL_D

Diameter Platform	Nomenclature	Connection type (Anti-rotational, Rotating or Angled)	Microdent Scanbody	Microdent Replica
 Ø4.10mm	1411 4.1 LAB (Laboratory)	E (Engaging - Anti-rotational)	CSUSBL41	
	1411 4.1 IO C (Short Intraoral)	NE (Non-engaging - Rotating) E (Engaging - Angled) NE (Non-engaging - Angled)	CSUSBI41C	CSURID41
	1411 IO L (Long Intraoral)		CSUSBI41L	
 Ø5.10mm	1412 5.1 LAB (Laboratory)	E (Engaging - Anti-rotational)	CSUSBL51	
	1412 5.1 IO C (Short Intraoral)	NE (Non-engaging - Rotating) E (Engaging - Angled) NE (Non-engaging - Angled)	CSUSBI51C	CSURID51
	1412 5.1 IO L (Long Intraoral)		CSUSBI51L	





FLOW DIRECT IMPLANT




PROCEDURE WITH DOME TI BASE


GALLERY: MICRODENT-1410_UNIVERSAL_TB

Diameter Platform	Nomenclature	Connection type (Anti-rotational, Rotating or Angled)	Ti Base Straight dome	Ti Base Dynamic dome
 Ø4.10mm	1411 3.5 LAB (Laboratory)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating)	Short CSUMDOM41H Short CSUMDOM41R	
	1411 3.5 IO C (Short Intraoral)		Long CSUMDOM41HL Long CSUMDOM41RL	
	1411 3.5 IO L (Long Intraoral)	E (Engaging - Angled) NE (Non-engaging - Angled)		Short CSUMDOMD41H Short CSUMDOMD41R Long CSUMDOMD41HL Long CSUMDOMD41RL
 Ø5.10mm	1412 4.2 LAB (Laboratory)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating)	Short CSUMDOM51H Short CSUMDOM51R	
	1412 4.2 IO C (Short Intraoral)		Long CSUMDOM51HL Long CSUMDOM51RL	
	1412 4.2 IO L (Long Intraoral)	E (Engaging - Angled) NE (Non-engaging - Angled)		Short CSUMDOMD51H Short CSUMDOMD51R Long CSUMDOMD51HL Long CSUMDOMD51RL

FLOW WITH STRAIGHT DOME TI BASE



FLOW WITH DYNAMIC DOME TI BASE



PROCEDURE


DIRECT MINI HEXAGONAL CAPITAL ABUTMENT

GALLERY: [MICRODENT-1430_Capitel_E_Mini_D](#)

Diameter Platform	Nomenclature	Connection type (Anti-rotational or Rotating)	Scanbody Mini Capitel	Replica Mini Capitel
Ø4.10mm	1431 Mini 4X LAB (Laboratory)		CCAPNSBL4XH	
	1431 Mini 4X IO C (Short Intraoral)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating)	CCAPNSBI4XHC	CCAPNRID4XH
	1431 Mini 4X IO L (Long Intraoral)		CCAPNSBI4XHL	
Ø5.10mm	1432 Mini 5X LAB (Laboratory)		CCAPNSBL51H	
	1432 Mini 5X IO C (Short Intraoral)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating)	CCAPNSBI51HC	CCAPNRID51H
	1432 Mini 5X IO L (Long Intraoral)		CCAPNSBI51HL	

FLOW

DIRECT MINI HEXAGONAL CAPITAL ABUTMENT





Microdent always recommends the use of the intraoral scanbody for a better fit between the prosthesis and the Mini Capitel Abutment using the one abutment one time technique.

PROCEDURE

TI BASE MINI HEXAGONAL CAPITAL ABUTMENT

GALLERY: [MICRODENT-1430_Capitel_E_Mini_TB](#)


Diameter Platform	Nomenclature	Connection type (Anti-rotational or Rotating)	Ti Base Straight dome
 Ø4.10mm	1431 Mini 4X LAB (Laboratory)		
	1431 Mini 4X IO C (Short Intraoral)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating)	CUTSNPC4XH
	1431 Mini 4X IO L (Long Intraoral)		
 Ø5.10mm	1432 Mini 5X LAB (Laboratory)		
	1432 Mini 5X IO C (Short Intraoral)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating)	CUTSNPC51H
	1432 Mini 5X IO L (Long Intraoral)		




FLOW TI BASE MINI HEXAGONAL CAPITAL ABUTMENT

Platform Diameter


- Ø4.20mm
- Ø5.10mm



Scanbodies on Mini Capital Transepithelial Abutment





Dome Ti Base for Mini Capital




PROCEDURE

DIRECT ROTATING MINI CAPITEL ABUTMENT

GALLERY: [MICRODENT-1440_Capitel_NE_Mini_D](#)

Diameter Platform	Nomenclature	Connection type (Rotating)	Scanbody Mini Capitel	Replica Mini Capitel
 Ø4.10mm	1441 Mini 4X LAB (Laboratory)		CCAPNSBL4X	
	1441 Mini 4X IO C (Short Intraoral)	NE (Non-engaging - Rotating)	CCAPNSBI4XC	CCAPNRID4X
	1441 Mini 4X IO L (Long Intraoral)		CCAPNSBI4XL	
 Ø5.10mm	1442 Mini 5X LAB (Laboratory)		CCAPNSBL51	
	1442 Mini 5X IO C (Short Intraoral)	NE (Non-engaging - Rotating)	CCAPNSBI51C	CCAPNRID51
	1442 Mini 5X IO L (Long Intraoral)		CCAPNSBI51L	

FLOW DIRECT ROTATING MINI CAPITEL ABUTMENT





Microdent always recommends the use of the intraoral scanbody for a better fit between the prosthesis and the Mini Capitel Abutment using the one abutment one time technique.

PROCEDURE

TI BASE ROTATING MINI CAPITAL ABUTMENT

GALLERY: [MICRODENT-1440_Capitel_NE_Mini_TB](#)

Diameter Platform	Nomenclature	Connection type (Rotating or Angled)	Ti Base Straight dome	Ti Base Dynamic dome
 Ø4.10mm	1441 Mini 4X LAB (Laboratory)	NE (Non-engaging - Rotating)	CUTSNPC4X	CUTSNDPC4X
	1441 Mini 4X IO C (Short Intraoral)			
	1441 Mini 4X IO L (Long Intraoral)			
 Ø5.10mm	1442 Mini 5X LAB (Laboratory)	NE (Non-engaging - Rotating)	CUTSNPC51	CUTSNDPC51
	1442 Mini 5X IO C (Short Intraoral)			
	1442 Mini 5X IO L (Long Intraoral)			




FLOW TI BASE ROTATING MINI CAPITAL ABUTMENT


Platform Diameter

■ Ø4.10mm

■ Ø5.10mm



Scanbodies on Mini Capitel Transepithelial Abutment



Mini Capitel Transepithelial Abutment





GENIUS CONICAL INTERNAL CONNECTION IMPLANT LIBRARY



PROCEDURE


DIRECT IMPLANT

GALLERY: MICRODENT-1450_GENIUS_D

Diameter Platform	Nomenclature	Connection type (Rotating or Angled)	Microdent Scanbody	Microdent Replica
 Ø3.00mm	1451 3.0 LAB (Laboratory)		CGN3SBL	
	1451 3.0 IO C (Short Intraoral)	NE (Non-engaging - Rotating) NE (Non-engaging - Angled)	CGN3SBIC	CGN3RID
	1451 3.0 IO L (Long Intraoral)		CGN3SBIL	
 Ø3.50mm Ø4.00mm Ø4.50mm Ø5.00mm	1451 3.5 LAB (Laboratory)		CGNSBL	
	1451 3.5 IO C (Short Intraoral)	NE (Non-engaging - Rotating) NE (Non-engaging - Angled)	CGNSBIC	CGNRID
	1451 3.5 IO L (Long Intraoral)		CGNSBIL	



FLOW DIRECT IMPLANT



PROCEDURE WITH DOME TI BASE

GALLERY: MICRODENT-1450_GENIUS_TB


Diameter Platform	Nomenclature	Connection type (Anti-rotational, Rotating or Angled)	Ti Base Straight dome	Ti Base Dynamic dome	
 Ø3.00mm		E (Engaging - Anti-rotational)	Short CGN3DOM3X00HC Short CGN3DOM3X02HC Short CGN3DOM3X03HC Short CGN3DOM3X04HC Short CGN3DOM3X05HC	Short CGN3DOMD3X00HC Short CGN3DOMD3X02HC Short CGN3DOMD3X03HC Short CGN3DOMD3X04HC Short CGN3DOMD3X05HC	
		E (Engaging - Angled)			
		1451 3.0 LAB (Laboratory)	NE (Non-engaging - Rotating)	Short CGN3DOM3X00RC Short CGN3DOM3X02RC Short CGN3DOM3X03RC Short CGN3DOM3X04RC Short CGN3DOM3X05RC	Short CGN3DOMD3X00RC Short CGN3DOMD3X02RC Short CGN3DOMD3X03RC Short CGN3DOMD3X04RC Short CGN3DOMD3X05RC
		1451 3.0 IO C (Short Intraoral)	NE (Non-engaging - Angled)		
		E (Engaging - Anti-rotational)	Long CGN3DOM3X00HL Long CGN3DOM3X02HL Long CGN3DOM3X03HL Long CGN3DOM3X04HL Long CGN3DOM3X05HL	Long CGN3DOMD3X00HL Long CGN3DOMD3X02HL Long CGN3DOMD3X03HL Long CGN3DOMD3X04HL Long CGN3DOMD3X05HL	
		E (Engaging - Angled)			
		1451 3.0 IO L (Short Intraoral)	NE (Non-engaging - Rotating)	Long CGN3DOM3X00RL Long CGN3DOM3X02RL Long CGN3DOM3X03RL Long CGN3DOM3X04RL Long CGN3DOM3X05RL	Long CGN3DOMD3X00RL Long CGN3DOMD3X02RL Long CGN3DOMD3X03RL Long CGN3DOMD3X04RL Long CGN3DOMD3X05RL
		NE (Non-engaging - Angled)			

PROCEDURE WITH DOME TI BASE


Diameter Platform	Nomenclature	Connection type (Anti-rotational, Rotating or Angled)	Ti Base Straight dome	Ti Base Dynamic dome
<div style="display: flex; flex-direction: column; align-items: center;"> <div style="width: 10px; height: 10px; background-color: #8bc34a; margin-bottom: 5px;"></div> <div style="text-align: center;"> <p>Ø3.50mm</p> <p>Ø4.00mm</p> <p>Ø4.50mm</p> <p>Ø5.00mm</p> </div> </div>		E (Engaging - Anti-rotational)	Short CGNDOM4X00HC Short CGNDOM4X02HC Short CGNDOM4X03HC Short CGNDOM4X04HC Short CGNDOM4X05HC	
		E (Engaging - Angled)		Short CGNDOMD4X00HC Short CGNDOMD4X02HC Short CGNDOMD4X03HC Short CGNDOMD4X04HC Short CGNDOMD4X05HC
	1452 3.5 LAB (Laboratory)	NE (Non-engaging - Rotating)	Short CGNDOM4X00RC Short CGNDOM4X02RC Short CGNDOM4X03RC Short CGNDOM4X04RC Short CGNDOM4X05RC	
	1452 3.5 IO C (Short Intraoral)	NE (Non-engaging - Angled)		Short CGNDOMD4X00RC Short CGNDOMD4X02RC Short CGNDOMD4X03RC Short CGNDOMD4X04RC Short CGNDOMD4X05RC
		E (Engaging - Anti-rotational)	Long CGNDOM4X00HL Long CGNDOM4X02HL Long CGNDOM4X03HL Long CGNDOM4X04HL Long CGNDOM4X05HL	
		E (Engaging - Angled)		Long CGNDOMD4X00HL Long CGNDOMD4X02HL Long CGNDOMD4X03HL Long CGNDOMD4X04HL Long CGNDOMD4X05HL
	1452 3.5 IO L (Short Intraoral)			
		NE (Non-engaging - Rotating)	Long CGNDOM4X00RL Long CGNDOM4X02RL Long CGNDOM4X03RL Long CGNDOM4X04RL Long CGNDOM4X05RL	
		NE (Non-engaging - Angled)		Long CGNDOMD4X00RL Long CGNDOMD4X02RL Long CGNDOMD4X03RL Long CGNDOMD4X04RL Long CGNDOMD4X05RL



FLOW WITH STRAIGHT DOME TI BASE





FLOW WITH DYNAMIC DOME TI BASE




PROCEDURE WITH DIRECT HEXAGONAL MINI CAPITEL

GALLERY: [MICRODENT-1430_Capitel_E_Mini_D](#)

Diameter Platform	Nomenclature	Connection type (Anti-rotational or Rotating)	Scanbody Mini Capitel	Replica Mini Capitel
 Ø3.00mm	1431 Mini 4X LAB (Laboratory)		CCAPNSBL4XH	
 Ø3.50mm Ø4.00mm Ø4.50mm Ø5.00mm	1431 Mini 4X IO C (Short Intraoral)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating)	CCAPNSBI4XHC	CCAPNRID4XH
	1431 Mini 4X IO L (Long Intraoral)		CCAPNSBI4XHL	





FLOW WITH DIRECT HEXAGONAL MINI CAPITEL



Microdent always recommends the use of the intraoral scanbody for a better fit between the prosthesis and the Mini Capitel Abutment using the one abutment one time technique.

PROCEDURE WITH TI BASE HEXAGONAL MINI CAPITEL

GALLERY: [MICRODENT-1430_Capitel_E_Mini_TB](#)

Diameter Platform	Nomenclature	Connection type (Anti-rotational or Rotating)	Ti Base Straight dome
 Ø3.00mm	1431 Mini 4X LAB (Laboratory)		CUTSNPC4XH
 Ø3.50mm Ø4.00mm Ø4.50mm Ø5.00mm	1431 Mini 4X IO C (Short Intraoral)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating)	
	1431 Mini 4X IO L (Long Intraoral)		CUTSNPC51H



FLOW WITH TI BASE HEXAGONAL MINI CAPITEL

Platform Diameter

- Ø3.00mm
- Ø3.50mm
- Ø4.00mm
- Ø4.50mm
- Ø5.00mm



Straight (heights)

- 1mm → CGN3CAPN4801H
→ CGNCAPN4801H
- 2mm → CGN3CAPN4802H
→ CGNCAPN4802H
- 3mm → CGN3CAPN4803H
→ CGNCAPN4803H
- 4mm → CGN3CAPN4804H
→ CGNCAPN4804H
- 5mm → CGN3CAPN4805H
→ CGNCAPN4805H



Scanbodies on Mini Capitel Transepithelial Abutment

Intraoral scanbody

Laboratory scanbody

Short

→ CCAPNSBI4XHC



→ CCAPNSBL4XH



Long

→ CCAPNSBI4XHL



Dome Ti Base for Mini Capitel



Straight


→ CUTSNPC4XH




PROCEDURE WITH DIRECT ROTATING MINI CAPITEL

GALLERY: [MICRODENT-1440_Capitel_NE_Mini_D](#)

Diameter Platform	Nomenclature	Connection type (Rotating)	Scanbody Mini Capitel	Replica Mini Capitel
 Ø3.00mm	1441 Mini 4X LAB (Laboratory)		CCAPNSBL4X	
 Ø3.50mm Ø4.00mm Ø4.50mm Ø5.00mm	1441 Mini 4X IO C (Short Intraoral)	NE (Non-engaging - Rotating)	CCAPNSBI4XC	CCAPNRID4X
	1441 Mini 4X IO L (Long Intraoral)		CCAPNSBI4XL	



FLOW WITH DIRECT ROTATING MINI CAPITEL





Microdent always recommends the use of the intraoral scanbody for a better fit between the prosthesis and the Mini Capitel Abutment using the one abutment one time technique.



PROCEDURE WITH TI BASE ROTATING MINI CAPITEL

GALLERY: [MICRODENT-1440_Capitel_NE_Mini_TB](#)


Diameter Platform	Nomenclature	Connection type (Rotating or Angled)	Ti Base Straight dome	Ti Base Dynamic dome
 Ø3.00mm	1441 Mini 4X LAB (Laboratory)	NE (Non-engaging - Rotating)	CUTSNPC4X	
 Ø3.50mm Ø4.00mm Ø4.50mm Ø5.00mm	1441 Mini 4X IO C (Short Intraoral)			
	1441 Mini 4X IO L (Long Intraoral)	NE (Non-engaging - Angled)		CUTSNDPC4X




FLOW WITH DOME TI BASE ON MINI CAPITEL

Platform Diameter


- Ø3.00mm
- Ø3.50mm
- Ø4.00mm
- Ø4.50mm
- Ø5.00mm



Scanbodies on Mini Capitel Transepithelial Abutment



Dome Ti Base for Mini Capitel



EK HEXAGONAL CONNECTION INTERNAL IMPLANT LIBRARY




PROCEDURE DIRECT IMPLANT

GALLERY: MICRODENT-1420_EKTOS_D




Diameter Platform	Nomenclature	Connection type (Anti-rotational, Rotating or Angled)	Microdent Scanbody	Microdent Replica
 Ø2.90mm	1421 2.9 LAB (Laboratory)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating) E (Engaging - Angled) NE (Non-engaging - Angled)	CEKSBL29	CEKRID29
	1421 2.9 IO C (Short Intraoral)		CEKSBI29C	
	1421 2.9 IO L (Long Intraoral)		CEKSBI29L	
 Ø3.50mm	1422 3.5 LAB (Laboratory)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating) E (Engaging - Angled) NE (Non-engaging - Angled)	CEKSBL35	CEKRID35
	1422 3.5 IO C (Short Intraoral)		CEKSBL35C	
	1422 3.5 IO L (Long Intraoral)		CEKSBL35L	
 Ø4.50mm	1423 4.5 LAB (Laboratory)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating) E (Engaging - Angled) NE (Non-engaging - Angled)	CEKSBL45	CEKRID45
	1423 4.5 IO C (Short Intraoral)		CEKSBI45C	
	1423 4.5 IO L (Long Intraoral)		CEKSBI45L	
 Ø5.70mm	1424 5.7 LAB (Laboratory)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating) E (Engaging - Angled) NE (Non-engaging - Angled)	CEKSBL57	CEKRID57
	1424 5.7 IO C (Short Intraoral)		CEKSBI57C	
	1424 5.7 IO L (Long Intraoral)		CEKSBI57L	

FLOW DIRECT IMPLANT




PROCEDURE WITH DOME TI BASE

GALLERY: MICRODENT-1450_EKTOS_TB


Diameter Platform	Nomenclature	Connection type (Anti-rotational, Rotating or Angled)	Ti Base Straight dome	Ti Base Dynamic dome
 Ø2.90mm	1421 2.9 LAB (Laboratory)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating) E (Engaging - Anti-rotational) NE (Non-engaging - Rotating)	Short CEKDOM29HC Short CEKDOM2901HC Short CEKDOM29RC Short CEKDOM2901RC Long CEKDOM29HL Long CEKDOM2901HL Long CEKDOM29RL Long CEKDOM2901RL	Short CEKDMD29HC Short CEKDMD2901HC Short CEKDMD29RC Short CEKDMD2901RC Long CEKDMD29HL Long CEKDMD2901HL Long CEKDMD29RL Long CEKDMD2901RL
	1421 2.9 IO C (Short Intraoral)			
	1421 2.9 IO L (Long Intraoral)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating) E (Engaging - Angled) NE (Non-engaging - Angled)		
 Ø3.50mm	1422 3.5 LAB (Laboratory)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating) E (Engaging - Anti-rotational) NE (Non-engaging - Rotating)	Short CEKDOM35HC Short CEKDOM3501HC Short CEKDOM35RC Short CEKDOM3501RC Long CEKDOM35HL Long CEKDOM3501HL Long CEKDOM35RL Long CEKDOM3501RL	Short CEKDMD35HC Short CEKDMD3501HC Short CEKDMD35RC Short CEKDMD3501RC Long CEKDMD35HL Long CEKDMD3501HL Long CEKDMD35RL Long CEKDMD3501RL
	1422 3.5 IO C (Short Intraoral)			
	1422 3.5 IO L (Long Intraoral)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating) E (Engaging - Angled) NE (Non-engaging - Angled)		
 Ø4.50mm	1423 3.5 LAB (Laboratory)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating) E (Engaging - Anti-rotational) NE (Non-engaging - Rotating)	Short CEKDOM45HC Short CEKDOM4501HC Short CEKDOM45RC Short CEKDOM4501RC Long CEKDOM45HL Long CEKDOM4501HL Long CEKDOM45RL Long CEKDOM34501RL	Short CEKDMD45HC Short CEKDMD4501HC Short CEKDMD45RC Short CEKDMD4501RC Long CEKDMD45HL Long CEKDMD4501HL Long CEKDMD45RL Long CEKDMD4501RL
	1423 3.5 IO C (Short Intraoral)			
	1423 3.5 IO L (Long Intraoral)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating) E (Engaging - Angled) NE (Non-engaging - Angled)		

PROCEDURE WITH DOME TI BASE


Diameter Platform	Nomenclature	Connection type (Anti-rotational, Rotating or Angled)	Ti Base Straight dome	Ti Base Dynamic dome
 Ø5.70mm	1424 3.5 LAB (Laboratory)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating) E (Engaging - Anti-rotational) NE (Non-engaging - Rotating)	Short CEKDOM57HC Short CEKDOM5701HC Short CEKDOM57RC Short CEKDOM5701RC Long CEKDOM57HL Long CEKDOM5701HL Long CEKDOM57RL Long CEKDOM5701RL	
	1424 3.5 IO C (Short Intraoral)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating) E (Engaging - Angled) NE (Non-engaging - Angled)		Short CEKDMD57HC Short CEKDMD5701HC Short CEKDMD57RC Short CEKDMD5701RC Long CEKDMD57HL Long CEKDMD5701HL Long CEKDMD57RL Long CEKDMD5701RL
	1424 3.5 IO L (Long Intraoral)			



FLOW WITH STRAIGHT DOME TI BASE







FLOW WITH DYNAMIC DOME TI BASE




PROCEDURE WITH DIRECT HEXAGONAL MINI CAPITEL

GALLERY: MICRODENT-1430_Capitel_E_Mini_D

Diameter Platform	Nomenclature	Connection type (Anti-rotational or Rotating)	Scanbody Mini Capitel	Replica Mini Capitel
 Ø2.90mm  Ø3.50mm  Ø4.50mm	1431 Mini 4X LAB (Laboratory)		CCAPNSBL4XH	
	1431 Mini 4X IO C (Short Intraoral)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating)	CCAPNSBI4XHC	CCAPNRID4XH
	1431 Mini 4X IO L (Long Intraoral)		CCAPNSBI4XHL	
 Ø5.70mm	1432 Mini 5X LAB (Laboratory)		CCAPNSBL51H	
	1432 Mini 5X IO C (Short Intraoral)	E (Engaging - Anti-rotational) NE (Non-engaging - Rotating)	CCAPNSBI51HC	CCAPNRID51H
	1432 Mini 5X IO L (Long Intraoral)		CCAPNSBI51HL	





FLOW WITH DIRECT HEXAGONAL MINI CAPITEL



Microdent always recommends the use of the intraoral scanbody for a better fit between the prosthesis and the Mini Capitel Abutment using the one abutment one time technique.

PROCEDURE WITH TI BASE HEXAGONAL MINI CAPITEL

GALLERY: [MICRODENT-1430_Capitel_E_Mini_TB](#)

Diameter Platform	Nomenclature	Connection type (Anti-rotational or Rotating)	Ti Base Straight dome
 Ø2.90mm  Ø3.50mm  Ø4.50mm	1431 Mini 4X LAB (Laboratory) 1431 Mini 4X IO C (Short Intraoral) 1431 Mini 4X IO L (Long Intraoral)	 E (Engaging - Anti-rotational) NE (Non-engaging - Rotating)	 CUTSNPC4X
 Ø5.70mm	1432 Mini 5X LAB (Laboratory) 1432 Mini 5X IO C (Short Intraoral) 1432 Mini 5X IO L (Long Intraoral)	 E (Engaging - Anti-rotational) NE (Non-engaging - Rotating)	 CUTSNPC51

FLOW WITH TI BASE HEXAGONAL MINI CAPITEL

Platform Diameter

- Ø2.90mm
- Ø3.50mm
- Ø4.50mm
- Ø5.70mm



Straight (heights)

- | | | |
|-----|---|--------------|
| 1mm | → | CEKCAPN2901H |
| | → | CEKCAPN3501H |
| | → | CEKCAPN4501H |
| | → | CEKCAPN5701H |
| 2mm | → | CEKCAPN2902H |
| | → | CEKCAPN3502H |
| | → | CEKCAPN4502H |
| | → | CEKCAPN5702H |
| 3mm | → | CEKCAPN2903H |
| | → | CEKCAPN3503H |
| | → | CEKCAPN4503H |
| 4mm | → | CEKCAPN2904H |
| | → | CEKCAPN3504H |
| | → | CEKCAPN4504H |



Scanbodies on Mini Capitel Transepithelial Abutment

Intraoral scanbody

Short

- CCAPNSBI4XHC
- CCAPNSBI4XHC
- CCAPNSBI51HC



Long

- CCAPNSBI4XHL
- CCAPNSBI4XHL
- CCAPNSBI51HL



Laboratory scanbody

- CCAPNSBL4XH
- CCAPNSBL4XH
- CCAPNSBL51H



Dome Ti Base for Mini Capitel





Straight

- CUTSNPC4XH
- CUTSNPC4XH
- CUTSNPC51H




PROCEDURE WITH DIRECT ROTATING MINI CAPITEL

GALLERY: MICRODENT-1440_Capitel_NE_Mini_D

Diameter Platform	Nomenclature	Connection type (Rotating)	Scanbody Mini Capitel	Replica Mini Capitel
 Ø2.90mm  Ø3.50mm  Ø4.50mm	1441 Mini 4X LAB (Laboratory)	NE (Non-engaging - Rotating)	CCAPNSBL4X	CCAPNRID4X
	1441 Mini 4X IO C (Short Intraoral)		CCAPNSBI4XC	
	1441 Mini 4X IO L (Long Intraoral)		CCAPNSBI4XL	
 Ø5.70mm	1442 Mini 5X LAB (Laboratory)	NE (Non-engaging - Rotating)	CCAPNSBL51	CCAPNRID51
	1442 Mini 5X IO C (Short Intraoral)		CCAPNSBI51C	
	1442 Mini 5X IO L (Long Intraoral)		CCAPNSBI51L	





FLOW WITH DIRECT ROTATING MINI CAPITEL



Microdent always recommends the use of the intraoral scanbody for a better fit between the prosthesis and the Mini Capitel Abutment using the one abutment one time technique.


PROCEDURE WITH TI BASE ROTATING MINI CAPITEL

GALLERY: [MICRODENT-1440_Capitel_NE_Mini_TB](#)

Diameter Platform	Nomenclature	Connection type (Rotating or Angled)	Ti Base Straight dome	Ti Base Dynamic dome
 Ø2.90mm  Ø3.50mm  Ø4.50mm	1441 Mini 4X LAB (Laboratory)	NE (Non-engaging - Rotating)	CUTSNPC4X	
	1441 Mini 4X IO C (Short Intraoral)			
	1441 Mini 4X IO L (Long Intraoral)	NE (Non-engaging - Angled)		CUTSNDPC4X
 Ø5.70mm	1442 Mini 4X LAB (Laboratory)	NE (Non-engaging - Rotating)	CUTSNPC51	
	1442 Mini 4X IO C (Short Intraoral)			
	1442 Mini 4X IO L (Long Intraoral)	NE (Non-engaging - Angled)		CUTSNDPC51



FLOW WITH TI BASE ROTATING MINI CAPITEL




SUMMARY

GALLERIES AND PROCEDURES

SYSTEM	GALLERY	PROCEDURE	PAGE
Microdent System	MICRODENT-1400_HEX_EXT_D	Direct implant	15
	MICRODENT-1400_HEX_EXT_TB	Direct Ti Base	17
	MICRODENT-1430_Capitel_E_Mini_D	Direct hexagonal transepithelial	20
	MICRODENT-1430_Capitel_E_Mini_TB	Ti Base hexagonal transepithelial	22
	MICRODENT-1440_Capitel_NE_Mini_D	Direct rotating transepithelial	24
	MICRODENT-1440_Capitel_NE_Mini_TB	Ti Base rotating transepithelial	26
Microdent Universal	MICRODENT-1410_UNIVERSAL_D	Direct implant	29
	MICRODENT-1410_UNIVERSAL_TB	Direct Ti Base	31
	MICRODENT-1430_Capitel_E_Mini_D	Direct hexagonal transepithelial	34
	MICRODENT-1430_Capitel_E_Mini_TB	Ti Base hexagonal transepithelial	36
	MICRODENT-1440_Capitel_NE_Mini_D	Direct rotating transepithelial	38
	MICRODENT-1440_Capitel_NE_Mini_TB	Ti Base rotating transepithelial	40
Microdent Genius	MICRODENT-1450_GENIUS_D	Direct implant	43
	MICRODENT-1450_GENIUS_TB	Direct Ti Base	45
	MICRODENT-1430_Capitel_E_Mini_D	Direct hexagonal transepithelial	49
	MICRODENT-1430_Capitel_E_Mini_TB	Ti Base hexagonal transepithelial	51
	MICRODENT-1440_Capitel_NE_Mini_D	Direct rotating transepithelial	53
	MICRODENT-1440_Capitel_NE_Mini_TB	Ti Base rotating transepithelial	55
Microdent EK	MICRODENT-1420_EKTOS_D	Direct implant	58
	MICRODENT-1420_EKTOS_TB	Direct Ti Base	60
	MICRODENT-1430_Capitel_E_Mini_D	Direct hexagonal transepithelial	64
	MICRODENT-1430_Capitel_E_Mini_TB	Ti Base hexagonal transepithelial	66
	MICRODENT-1440_Capitel_NE_Mini_D	Direct rotating transepithelial	68
	MICRODENT-1440_Capitel_NE_Mini_TB	Ti Base rotating transepithelial	70





CUSTOMER SERVICE

Telephone 93 844 76 50

Hours: Monday to Thursday from 9am to 6pm.

Friday from 9am to 3pm

PLACING ORDERS

Telephone 93 844 76 50

E-mail: implant@microdentsystem.com

www.microdentsystem.com

Sales representative for your area.

DELIVERY TIMES

Delivery Mainland Spain before 7pm the next day for orders placed before 2pm (Iberian Peninsula).

Delivery between 24-48h for orders placed before 2pm (Canary Islands).

Consult same day delivery options (available only in Barcelona province).

DISPATCH OF PRODUCTS

Always for orders placed before 2pm.

Tell us your chosen time when placing your order:

- 8.30am Service.
- Preferential: delivery from 8.30am to 10am
- Express: delivery from 10am to 1pm.
- 7pm service: delivery before 7pm. (Default times).

*Available for the whole of mainland Spain. Consult for the Canary Islands.

RETURNS

Returns will not be accepted after a period of 30 days from the delivery of the goods.

All material returns must be accompanied by the product return form duly completed and enclosing a copy of the delivery note.

INVOICING AND METHOD OF PAYMENT

Bank transfer within 30 days of the invoice date.

Consult other payment methods and facilities with our customer service department or the sales representative for your area.

* Microdent reserves the right to make any changes without prior notice.

Important:

- The use of attachments not supplied by Microdent can compromise the stability and fixation of the prosthesis as well as cause loosening or fracture. It also renders the warranty for our implants null and void.

- Consult availability of references according to country as there may be products not marketed in your area.



