

For the correct use of our products, it is recommended that professionals receive training on their specific applications. Implant Microdent System S.L.U. regularly organizes training courses that provide information on the characteristics of each prosthetic component and CAD/CAM support solution, as well as their correct technical application. This information is also available on the company website, www.microdentsystem.com, and in our various technical publications.

INTENDED USE

Implant Microdent System S.L.U. offers the products necessary for the application and use of its prosthetic rehabilitation systems and digital solutions.

PRODUCT DESCRIPTION

For proper product selection, we recommend consulting our published technical catalogs. The prosthetic components and digital support offered by Implant Microdent System S.L.U. are manufactured from stainless steel, high-strength titanium, and high-strength polymer (POM).

PRODUCT PRESENTATION

Packaging

Implant Microdent System S.L.U. guarantees that all its products undergo a controlled manufacturing, inspection, and cleaning process before packaging. The products are available in different packaging options, each properly identified with its corresponding technical label.

Product supply state







Microdent prosthetic products are supplied in a clean, but not sterile, condition. If the technical workflow requires disinfection or deep cleaning before use in the prosthetic process, it is recommended to follow standard laboratory protocols.

If a thermal cleaning process is required due to workflow requirements, the specific characteristics of the equipment used must be taken into account. As a technical reference for thermoplastic and metallic materials, the thermal resistance parameters are:

- Maximum temperature: 134°C
- Exposure time: 3 minutes.

LABELLING SYMBOLS

Product label shows the following indications:

	Reference Number.
	Please follow the instructions for use.
	Lot number.
	Precaution, consult the warnings.
	Manufacturer.
	Unique identifier code.

WARNINGS AND PRECAUTIONS

Microdent products are part of a comprehensive system designed for compatibility with Implant Microdent System S.L.U. systems. The use of non-Microdent components may result in suboptimal performance.

These products must be used exclusively for their intended purpose. Operators or users handling these products are advised of the risk of sensitivity to the materials described in the "PRODUCT DESCRIPTION" section.

Before use, it is the responsibility of the technician to check:

- The condition of the packaging and that the reference number matches the product required for the prosthetic work.
- The integrity of the product to verify that it is in optimal condition.
- Maintenance and cleaning:
 - For proper product storage, the following guidelines must be followed:
 - Do not use chlorine-based or acid-based cleaning or disinfecting agents.
 - Avoid products containing aldehydes.
 - Thoroughly remove any residue from the surface or interior of the product using a nylon bristle brush.
 - For products with edges or adjustment areas, avoid impacts or contact with other elements that could damage the surface or cutting edges.
- After cleaning, rinse with distilled water and dry thoroughly. Improper drying may cause surface oxidation spots.

Non-medical products are subject to wear and tear from use. It is the user's responsibility to periodically replace them to ensure the accuracy of prosthetic adjustments. Implant Microdent System S.L.U. recommends always having replacement units available.

LIABILITY, SAFETY AND WARRANTY

If the packaging or product is damaged during transport, Implant Microdent System S.L.U. will replace it free of charge. Evidence of improper handling of the product exempts the company from the aforementioned liability.

The warranty applies to Implant Microdent System S.L.U. products provided they have been used according to the instructions for use and technical recommendations provided. To return the products mentioned in these instructions, the procedures established in our terms and conditions of sale and delivery must be followed.

Any technical anomaly must be carefully analyzed to determine its origin and any deviations in the work result, and must be properly documented.

Manufacturer by: Implant Microdent System S.L.U.
Pol. Ind. Can Magre C/. Carles Buigas, 1
08187 Sta. Eulàlia de Ronçana (Barcelona)